

# Empore™

## SDB-RPS Reverse Phase Sulfonated Extraction Disks for Environmental Analysis

*For use with 47 mm and 90 mm extraction apparatus*

### Instructions for Use

Note: Empore Sample Preparation Products are intended for solid phase extraction during scientific research only. These products are not intended for use in medical devices or in assessment and treatment of clinical patients.

### General Product Characteristics

#### **Description:**

Empore SDB-RPS (styrene-divinylbenzene-reverse phase sulfonate) Disks are a unique product for the solid phase extraction of polar and nonpolar analytes from aqueous samples. A proprietary process is used to entrap adsorbent particles into a matrix of inert PTFE to create a mechanically stable sorbent disk.

#### **Formulation:**

90% or greater adsorbent particle  
10% or less PTFE

#### **Product Characteristics:**

*Thickness:* 0.5 mm ± 0.05 mm  
*SPE Flow Rate:* < 10 min/L DI H<sub>2</sub>O @ 25°C  
@ 20 in. Hg (47 mm disk)  
*Solvents:* Compatible with all organic solvents  
*pH:* Stable between 1 and 14.

# Suggested Application Procedures

## General: Water Analysis

Empore™ Extraction Disks provide an efficient alternative to liquid/liquid extraction for sample preparation. A proprietary process is used to entrap adsorbent particles into a matrix of inert PTFE to create a mechanically stable sorbent disk. The disks can be used for purification and concentration of analytes for analysis. Advantages of Empore Extraction Disks include rapid filtration, reduced solvent usage and a reduction of analytical interferences.

**The enclosed instructions are general guidelines for use. Sample volume, solvent type and conditioning may be changed to adapt to specific methods and analytes as needed.**


## Extraction Equipment

- 47 mm or 90 mm Empore Extraction Disks.
- 47 mm or 90 mm glass filtration apparatus.
- Vacuum source.

## Sample Preparation

- Microbiological growth can be retarded by lowering sample pH to 2.
- Filter Aid 400 (Cat. # FA400) and/or prefiltration may be helpful if the sample contains excessive suspended solids.

# Suggested Vacuum Apparatus



Note: Suggested solvent volumes will vary according to the disk diameter and the amount of Filter Aid 400 filter material. A general guide for solvent volumes is to completely cover the disk and bed of filter material, such that 2 - 3 mm of solvent is above the surface. Repeat with second aliquot.

## Extraction Disk Conditioning

Disk conditioning is critical for a successful extraction. Conditioning provides a good interface between the sorbent and the sample matrix. FAILURE TO CONDITION THE EXTRACTION DISKS PROPERLY WILL RESULT IN ERRATIC AND LOW RECOVERY.

1. Center the extraction disk on the base of the filtration apparatus and clamp the reservoir on the top of the disk.\*
2. Wash the disk with 10 mL of acetone and apply vacuum to dry the disk.
3. Wash the disk with 10 mL isopropanol. Apply vacuum and dry the disk.
4. Add 10 mL methanol to the disk. Apply vacuum and pull approximately 1 mL through the disk. Vent the vacuum and allow the disk to soak for 30 seconds.
5. Apply vacuum and draw methanol through the disk leaving a small amount of methanol on surface.
6. Add 10 mL reagent grade water to the reservoir and draw the water through the disk until the water surface just covers the disk surface.

Note: When using solvents or other chemicals, be sure to read and follow the manufacturer's precautions and directions for use.

If disk becomes dry while conditioning with methanol or water, repeat steps 4 through 6.

- \* Place a vial in the vacuum apparatus to collect and dispose of wash and conditioning solvents. Remove vial prior to sample extraction.

## Sample Extraction

- Pour the sample into the reservoir and apply vacuum. Flow rate is dependent on vacuum source and solids content of the sample. However, recoveries are not affected by flow rate.
- After sample extraction is complete, remove residual water from the disk by applying vacuum to the disk for approximately 5-20 minutes.

## Sample Elution

Two elutions with 10 mL solvent are recommended.

- Place tip of filter base into the collection vessel (see diagram).
- Add 10 mL elution solvent to sample container carefully rinsing the sides. Transfer solvent from sample container to reservoir with a pipet washing the walls of the reservoir in the process.
- Apply vacuum and pull approximately 1 mL elution solvent through the disk. Vent the vacuum and allow the disk to soak 30 seconds before reapplying vacuum to dry the disk.
- Repeat this process with a second aliquot of eluting solvent.
- The sulfonic acid groups may require elution solvents containing pH modifiers or buffers such as ammonium hydroxide to elute some cationic analytes.

# General Information

## Handling and Storage

The disks may be handled in the same manner as any filter membrane. Because of the adsorptive properties of the disk, desiccator storage away from laboratory air at room temperature is recommended.

## Recommended Usage

Empore™ Extraction Disks are used in a manner similar to membrane filters. Filtration equipment is available from a number of different suppliers and include in-line filter holders, glass filtration apparatus, and multiple filtration manifolds. Buchner funnels are not recommended.

## Applications

Empore™ SDB-RPS (styrene-divinylbenzene-reverse phase sulfonate) Extraction Disks are used for extraction of polar and nonpolar organic compounds from water samples or soil extracts.

## **Important Notice**

3M MAKES NO WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Since a variety of factors can affect the use and performance of a 3M product in a particular application, some of which are uniquely within the user's knowledge and control, user is responsible for determining whether or not the 3M product is fit for a particular purpose and suitable for user's method of application.

## **Limitation of Remedies and Liability**

If the 3M product is proved to be defective, THE EXCLUSIVE REMEDY, AT 3M'S OPTION, SHALL BE TO REFUND THE PURCHASE PRICE OR TO REPAIR OR REPLACE THE DEFECTIVE PRODUCT. 3M shall not otherwise be liable for loss or damages, whether direct, indirect, special, incidental or consequential, regardless of the legal theory asserted, including, but not limited to, contract, negligence, warranty or strict liability.

## **Warranty**

Seller warrants its equipment against defects in workmanship and material for a period of 12 months from date of shipment from the factory under normal use and service and otherwise when such equipment is used in accordance with instructions furnished by Seller and for purposes disclosed in writing at the time of purchase, if any. Any unauthorized alteration or modification of the equipment by Buyer will void this warranty. Seller's liability under this warranty shall be limited to the replacement or repair, F.O.B., point of manufacture, of any defective equipment or part which, having been returned to the factory, transportation charges prepaid, has been inspected and determined by Seller to be defective. THIS WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EITHER EXPRESSED OR IMPLIED, AS TO DESCRIPTION, QUALITY, MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE OR USE, OR ANY OTHER MATTER. Under no circumstances shall Seller be liable to Buyer or any third party for any loss of profits or other direct or indirect costs, expenses, losses or consequential damages arising out of or as a result of any defects in or failure of

its products or any part or parts thereof or arising out of or as a result of parts or components incorporated in Seller's equipment but not supplied by the Seller.

### **Technical Information and Data**

Technical information and data, recommendations, and other statements provided by 3M are based on information, tests, or experience which 3M believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. No license under any 3M or third party intellectual property rights is granted or implied with this information.

### **Safety and Disposal**

Always follow safe handling practices while processing Empore products as applicable to your use. Always wear appropriate personal protective equipment while processing Empore products as applicable to your use. Dispose of used Empore products per requirements of your process. Alteration or addition of other materials to Empore products may require unique disposal protocols. Always dispose of used product in accordance with federal, state and local requirements.

3M is a trademark of 3M Company.  
Empore is a trademark of 3M Company used under license.

©2009 3M. All Rights Reserved. Issued: 3/09

## **Contact Information**

**To place orders or reach  
Customer or Technical Service:**

Toll Free: +1 888 509 5330

3130 Lexington Ave. S.  
Eagan, MN 55121

34-7054-6369-2

