

IMG Mouse Microglial Cell Line

Immortalized Cell Line

Cat. # SCC134

Pack size: $\geq 1 \times 10^6$

viable cells/vial

Store in liquid nitrogen

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.
NOT FOR HUMAN OR ANIMAL CONSUMPTION.
THIS PRODUCT CONTAINS GENETICALLY MODIFIED ORGANISMS.



Certificate of Analysis

page 1 of 3

Background

Microglia are the primary immune cells of the central nervous system (CNS), and account for 10-15% of all cells found within the brain. Microglial cells play an important role in maintaining the health of the CNS by the removal of pathogens, infectious agents and damaged cells through phagocytosis. Microglial cells have also been extensively studied for their harmful roles in a variety of neurodegenerative disease including Alzheimer's disease, Parkinson's disease and amyotrophic lateral sclerosis.

Immortalized **MicroG**lial Cell Line or IMG is an immortalized microglial cell line derived from adult murine brain infected with the v-raf/v-myc retrovirus. IMG cells express markers that are specific to adult primary microglial cells (F4/80 and CD11b) as assessed by flow cytometry and are responsive to pro- and anti-inflammatory signals¹. They have been shown to phagocytose foreign particles and misfolded proteins such as the amyloid-beta peptide¹. IMG cells are rapidly dividing cells that are typically small and round with few or no processes when cultured in medium containing serum.

Storage and Handling

IMG Mouse Microglial Cell Line should be stored in liquid nitrogen. The cells can be cultured for at least 10 passages after initial thawing without significantly affecting the cell marker expression and functionality.

Quality Control Testing

- Each vial contains $\geq 1 \times 10^6$ viable cells.
- Cells are tested negative for infectious diseases by a Mouse Essential CLEAR panel by Charles River Animal Diagnostic Services.
- Cells are verified to be of mouse origin and negative for inter-species contamination from rat, chinese hamster, Golden Syrian hamster, human and non-human primate (NHP) as assessed by a Contamination Clear panel by Charles River Animal Diagnostic Services
- Cells are negative for mycoplasma contamination.

Data




Figure 1. IMG cells exhibit similar morphology to primary microglia (A) and express microglial markers F4/80 (B) and CD11b (C). IMG cells were stained with anti-F4/80 (mouse), APC clone BM8.1 (Cat No. MABF1523; filled trace) (B). Rat IgG2a negative control, clone 2A3, APC conjugate (Cat. No. MABF1774; open trace) was used as the isotype control. IMG cells were stained with anti-CD11b (human/mouse), APC clone M1/70 (Cat. No. MABF520; filled trace) (C). Rat IgG2b negative control, clone LTF-2, APC conjugate (Cat. No. MABF1772; open trace) was used as the isotype control.

References

1. McCarthy RC, Lu DY, Alkhateeb A, Gardeck AM, Lee CH and Wessling-Resnick M (2016). Characterization of a novel adult murine immortalized microglial cell line and its activation by amyloid-beta. *J. Neuroinflammation* 13: 21

SPECIES LEGEND: H Human Ca Canine M Mouse R Rat Rb Rabbit B Bovine P Porcine WR Most Common Vertebrates

Please visit www.millipore.com for additional product information and references.

Submit your published journal article, and earn credit toward future purchases. Visit www.millipore.com/publicationrewards to learn more!

Protocols

Thawing Cells

1. Do not thaw the cells until the recommended medium is on hand. Cells can grow on normal tissue cultureware surfaces without any additional coating.

IMG Expansion Medium: Cells are thawed and expanded in High Glucose DMEM (Sigma Cat. No. D6546), 10% FBS (Cat. No. ES-009-B), 1X L-Glutamine (Cat. No. TMS-002-C) and 1X Penicillin-Streptomycin Solution (Cat. No. TMS-AB2-C) (optional).

2. Remove the vial of frozen IMG cells from liquid nitrogen and incubate in a 37°C water bath. Closely monitor until the cells are completely thawed. Maximum cell viability is dependent on the rapid and complete thawing of frozen cells.

IMPORTANT: Do not vortex the cells.

3. As soon as the cells are completely thawed, disinfect the outside of the vial with 70% ethanol. Proceed immediately to the next step.
4. In a laminar flow hood, use a 1 or 2 mL pipette to transfer the cells to a sterile 15 mL conical tube. Be careful not to introduce any bubbles during the transfer process.
5. Using a 10 mL pipette, slowly add dropwise 9 mL of IMG Expansion Medium (Step 1 above) to the 15 mL conical tube.

IMPORTANT: Do not add the entire volume of media all at once to the cells. This may result in decreased cell viability due to osmotic shock.

6. Gently mix the cell suspension by slowly pipetting up and down twice. Be careful not to introduce any bubbles.

IMPORTANT: Do not vortex the cells.

7. Centrifuge the tube at 300 x g for 2-3 minutes to pellet the cells.
8. Decant as much of the supernatant as possible. Steps 5-8 are necessary to remove residual cryopreservative (DMSO).
9. Resuspend the cells in 10-15 mL of IMG Expansion Medium.
10. Transfer the cell mixture to a T75 tissue culture flask.
11. Incubate the cells at 37°C in a humidified incubator with 5% CO₂.
12. The next day, exchange the medium with 10-15 mL of fresh IMG Expansion Medium. Exchange with fresh medium every two to three days thereafter.
13. When the cells are approximately 80-85% confluent, they can be dissociated with Accutase (Cat. No. SCR005) or trypsin-EDTA (Cat. No. SM-2003-C) and further passaged or, alternatively, frozen for later use.

Subculturing Cells

1. Carefully remove the medium from the T75 tissue culture flask containing the confluent layer of IMG cells.
 2. Rinse the T75 flask twice with 10 mL 1X PBS. Aspirate after each rinse.
 3. Apply 3-5 mL of Accutase or trypsin-EDTA solution and incubate in a 37°C incubator for 3-5 minutes.
 4. Inspect the flask and ensure the complete detachment of cells by gently tapping the side of the flask with the palm of your hand.
 5. Add 8 mL of IMG Expansion Medium to the plate.
 6. Gently rotate the flask to mix the cell suspension. Transfer the dissociated cells to a 15 mL conical tube.
 7. Centrifuge the tube at 300 x g for 3-5 minutes to pellet the cells.
 8. Discard the supernatant, then loosen the cell pellet by tapping the tip of the tube with a finger.
 9. Apply 2 mL of IMG Expansion Medium to the conical tube and resuspend the cells thoroughly.
- IMPORTANT: Do not vortex the cells.**
10. Count the number of cells using a hemocytometer.
 11. Plate the cells to the desired density (typical split ratio is 1:5 – 1:7).

Cryopreservation of Cells

IMG Mouse Microglial Cell Line may be frozen in the expansion medium plus 8-10% DMSO using a Nalgene slow freeze Mr. Frosty container.

RESTRICTED USE AGREEMENT
(subject to local law)

THIS PRODUCT MAY ONLY BE USED FOR RESEARCH PURPOSES, WHICH IS FURTHER DEFINED BELOW. BY OPENING THIS PRODUCT, YOU (“PURCHASER”) HEREBY REPRESENT THAT YOU HAVE THE RIGHT AND AUTHORITY TO LEGALLY BIND YOURSELF AND/OR YOUR EMPLOYER, AS APPLICABLE, AND CONSENT TO BE LEGALLY BOUND BY THE TERMS OF THIS RESTRICTED USE AGREEMENT. IF YOU DO NOT AGREE TO COMPLY WITH THESE TERMS, YOU MAY NOT OPEN OR USE THE PRODUCT AND YOU MUST CALL MILLIPORESIGMA (“SELLER”) CUSTOMER SERVICE (1-800-645-5476) TO ARRANGE TO RETURN THE PRODUCT FOR A REFUND.

“Product” means IMG Mouse Microglial Cell Line (SCC134)

“Research Purposes” means any internal *in vitro* research use and specifically excludes the following uses of whatever kind or nature:

- Re-engineering or copying the Product
- Making derivatives, modifications, or functional equivalents of the Product
- Obtaining patents or other intellectual property rights claiming use of the Product
- Using the Product in the development, testing, or manufacture of a Commercial Product
- Using the Product as a component of a Commercial Product
- Reselling or licensing the Product
- Using the Product in clinical or therapeutic applications including producing materials for clinical trials
- Administering the Product to humans
- Using the Product in collaboration with a commercial or non-academic entity

“Commercial Product” means any product intended for: (i) current or future sale; (ii) use in a fee-for-service; or (iii) any diagnostic, clinical, or therapeutic use.

Access to the Product is limited solely to PURCHASER’s officers, employees, and students who need to use the Product for Research Purposes. PURCHASER shall comply with all applicable laws in its use and handling of the Product and shall keep it under reasonably safe and secure conditions to prevent unauthorized use or access.

These restrictions will remain in effect for as long as PURCHASER possesses the Product.

PLEASE CONTACT licensing@emdmillipore.com PRIOR TO PURCHASE FOR ANY USE OF THE PRODUCT OUTSIDE OF THIS RESTRICTED USE AGREEMENT.

GMO

This product contains genetically modified organisms.
Este producto contiene organismos genéticamente modificados.
Questo prodotto contiene degli organismi geneticamente modificati.
Dieses Produkt enthält genetisch modifizierte Organismen.
Ce produit contient des organismes génétiquement modifiés.
Dit product bevat genetisch gewijzigde organismen.
Tämä tuote sisältää geneettisesti muutettuja organismeja.
Denna produkt innehåller genetiskt ändrade organismer.

■ antibodies ■ Multiplex products ■ biotools ■ cell culture ■ enzymes ■ kits ■ proteins/peptides ■ siRNA/cDNA products

Please visit www.millipore.com for additional product information, test data and references

EMD Millipore Corporation, 28820 Single Oak Drive, Temecula, CA 92590, USA 1-800-437-7500

Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. Not for human or animal consumption. Purchase of this Product does not include any right to resell or transfer, either as a stand-alone product or as a component of another product. Any use of this Product for purposes other than research is strictly prohibited.

Millipore®, the M mark, Upstate®, Chemicon®, Linco® and all other trademarks, unless specifically identified above in the text as belonging to a third party, are owned by Merck KGaA, Darmstadt. Copyright ©2008-2017 Merck KGaA, Darmstadt. All rights reserved.



We Buy 100% Certified
Renewable Energy