## Instructions for Use

## Reticulocyte Stain

## Procedure No. R4132

## Intended Use

Reticulocyte Stain is intended for the identification of reticulocytes in human blood films. Reticulocyte Stain reagents are for "In Vitro Diagnostic Use." For professional use only. The data obtained from this manual qualitative procedure serves only as an aid to diagnosis of hemolytic anemia or to assess bone marrow function and should be reviewed in conjunction with other clinical diagnostic tests or information.

In 1949, Brecher<sup>1</sup> introduced the new methylene blue method for identification of reticulocytes based on precipitation of ribosomal RNA by the cationic dye. This has now replaced other methods and is the recognized procedure for the manual quantitation of reticulocytes in peripheral blood.<sup>2</sup>

Sigma-Aldrich provides a stable solution for the staining of reticulocytes in blood films. The blood and reticulocyte staining solution are mixed and incubated briefly at room temperature. Wedge or spun smears are made on microscope slides, air dried and evaluated under oil immersion on a light microscope. A red blood cell scoring positive would be observed containing two or more blue-stained granules.

## Reagents

**Reticulocyte Stain** (Cat. No. R4132-120ML)

New methylene blue, 1% (w/v), C.I. 52030, sodium chloride, 0.9% (w/v), and potassium oxalate, 1.6% (w/v)

## Special Materials Required but Not Provided

- Microscope
- Blood smearing instrument or Cytocentrifuge
- Glass tubes, 10x75 mm or 12x75 mm
- Laboratory droppers or Pasteur pipets
- Microscope slides and coverslips

## Storage and Stability

Store Reticulocyte Stain at room temperature (18–26°C). Reagent is stable until expiration date.

## Preparation

Reticulocyte Stain is supplied ready for use.

## Precautions

This IVD is intended for in vitro diagnostic use in a clinical laboratory environment. This IVD is for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

## Procedure

## Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Blood may be collected in standard clinical laboratory evacuated tubes. All routine anticoagulants are acceptable (e.g., heparin, citrate and oxalate). The tripotassium salt of ethylenediamine tetraacetic acid (K3 EDTA) is the anticoagulant of choice. If procedure is not carried out within 2–3 hours of collection, store blood at 4°C. Blood should be warmed to room temperature and mixed thoroughly prior to staining. Blood over 24-hours old is not recommended for use.

## Notes

It is recommended that blood smears prepared from healthy donors be processed along with patient samples as normal controls.

A small amount of precipitate may form in the Reticulocyte Stain. If precipitate is noticed, filter through laboratory grade filter paper.

## Procedure

1. Add three drops of thoroughly mixed blood and two drops of Reticulocyte Stain, to a glass tube and mix well. Let stand 10 minutes at room temperature (18–26°C).
2. Prepare a conventional wedge or spun smear and air dry at least 15 minutes. Counterstaining is not recommended.
3. Coverslip and examine microscopically.

## Performance Characteristics

Stained blood films are evaluated for the presence or absence of reticulocytes. A reticulocyte is considered any red blood cell containing two or more blue-stained particles.<sup>2</sup>

Using a 100x oil immersion objective and a 10x ocular, randomly pick areas of the film where red cells are close to each other but do not touch or overlap. Count 1000 red blood cells including reticulocytes. The proportion of reticulocytes may be calculated as:

$$\text{Reticulocyte Count (\%)} = \frac{\text{Total Number of Reticulocytes}}{10}$$

Normal adult values at sea level are 1.0 ± 0.5%.<sup>2</sup>

Due to variability in the hematocrit, it may be necessary to correct the observed reticulocyte count to a normal hematocrit of 45%:

$$\text{Corrected Reticulocyte Count (\%)} = \text{Observed Count} \times \frac{\text{Measured Hematocrit (\%)}}{45\%}$$

In the normal non-anemic patient, reticulocytes are found in circulating blood on the fourth day of maturation following three days of maturation within the marrow. Factors increasing erythropoiesis may shorten bone marrow maturation time while lengthening blood maturation time. This shift results in large numbers of "shift" reticulocytes circulating in the blood and these should not be considered when quantitating reticulocytes as a reflection of red blood cell production.<sup>2,3</sup> Shift cells may be detected by a Wright-stained film<sup>3</sup> and a reticulocyte production index<sup>2</sup> established on the corrected reticulocyte count and hematocrit:

$$\text{Reticulocyte Production Index} = \frac{\text{Corrected Reticulocyte Count (\%)}}{\text{Expected Maturation Time (Days)}}$$

With expected maturation time as follows:

Days	Hematocrit
0	45%
1.5	35%
2	25%
3	15%

A reticulocyte production index greater than or equal to 3 is considered normal, while an index of less than 2 is below normal.<sup>2</sup>

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

## Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
R4132	Reticulocyte Stain	Reticulocytes	3 of 3	3 of 3	3 of 3	3 of 3

## Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

## R4132:

Not a hazardous substance or mixture.

**If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.**

## Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer

## References

1. Brecher G: New methylene blue as a reticulocyte stain. Am J Clin Pathol 19:895, 1949
2. Koepke JA. Update on reticulocyte counting. Lab Med 1999;30:339–343
3. Perrotta AL, Finch CA: The poly chromatophilic erythrocyte. Am J Clin Pathol 57:471, 1972

## Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](http://SigmaAldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](http://SigmaAldrich.com/techservice).

## Revision History

Rev. 5.0	2016
Rev. 6.0	2022
Rev. 7.0	2022

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Updated literature references. Added a revision history table. Added Warnings and Hazards.

Sigma-Aldrich, Inc.,  
3050 Spruce Street,  
St. Louis, MO 63103 USA  
an affiliate of Merck KGaA,  
Darmstadt, Germany  
+1(314) 771-5765

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Merck KGaA,  
64271 Darmstadt,  
Germany

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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MDSS GmbH