


Product Information

Bryostatin 1

Product Number **B 7431**
 Storage Temperature $-20\text{ }^{\circ}\text{C}$

CAS# 83314-01-6



Product Description

Molecular Formula: $\text{C}_{47}\text{H}_{68}\text{O}_{17}$

Molecular Weight: 905.0

Supplied as white powder or colorless film

Purity: Approximately 99% by HPLC

Bryostatin 1 is a macrolactone isolated from the marine bryozoan, *Bugula neritina*, that modulates protein kinase C (PKC) activity. Bryostatins compete with the phorbol esters and diacylglycerols for binding to the cysteine-rich domains CRD1 and CRD2 of PKC. However, while the phorbol esters are tumor promoters, the bryostatins have antineoplastic activity. Bryostatin 1 induces an initial rapid activation and autophosphorylation of PKC that results in the translocation of the PKC enzyme to the membrane. PKC is subsequently down-regulated by ubiquitination and proteolysis.¹

The effects of bryostatin 1 are cell type and PKC isoform dependent. For example, in LNCaP human prostate cancer cells, bryostatin 1 induces a transient activation and membrane translocation of PKC α followed by down-regulation, while it induces a prolonged translocation of PKC δ and PKC ϵ to non-nuclear membranes.² Other studies have also indicated that PKC δ is not down-regulated by bryostatin 1.³

Through its activation of PKC, bryostatin-1 induces the phosphorylation of other kinases and regulatory proteins. For example, in NB4 promyelocytic leukemia cells bryostatin 1 induces cell differentiation and the hyperphosphorylation of an abnormal p53 via the MAP kinase signaling pathway. Hyperphosphorylated p53 is subsequently ubiquitinated and degraded.⁴ Bryostatin 1 also induces the phosphorylation of the oncogene Bcl-2 by PKC α . In the acute lymphoblastic leukemia cell line, Reh, the phosphorylation of Bcl-2 leads to its ubiquitin-mediated degradation.⁵ However, in mouse embryonic fibroblasts (MEF), the phosphorylation of Bcl-2 protects against heat-induced apoptosis.⁶

Bryostatin 1 sensitizes cancer cells to the cytotoxic effects of anti-cancer agents. For example, chronic preexposure of leukemia cell lines to bryostatin 1 down-regulates PKC and enhances the cytotoxicity of 1- β -D-arabinofuranosylcytosine (Ara-C).⁷ In Reh cells bryostatin 1 enhances the the cytotoxicity of auristatin PE and vincristine, perhaps via the down-regulation of Bcl-2.⁸

Preparation Instructions

Bryostatin 1 is soluble in DMSO and ethanol.

Storage/Stability

Store tightly sealed at $-20\text{ }^{\circ}\text{C}$. Bryostatin 1 binds to glass and plastic surfaces in aqueous solutions.

References

1. Mutter, R., and Wills, M., Chemistry and clinical biology of the bryostatins. *Bioorg. Med. Chem.*, **8**, 1841-1860 (2000).
2. Gschwend, J.E., et al., Bryostatin 1 induces prolonged activation of extracellular regulated protein kinases and apoptosis of LNCaP cancer cells overexpressing protein kinase c alpha. *Mol. Pharmacol.*, **57**, 1224-1234 (2000).
3. Lorenzo, P.I., et al., The catalytic domain of protein kinase C δ confers protection from down-regulation induced by bryostatin 1. *J. Biol. Chem.*, **272**, 33338-33343 (1997).
4. Song, X.D., et al., Mitogen-activated protein kinase is involved in the degradation of p53 protein in the Bryostatin-1-induced differentiation of the acute promyelocytic leukemia NB4 cell line. *J. Biol. Chem.*, **274**, 1677-1682 (1999).
5. Wall, N.R., et al., Bryostatin 1 induces ubiquitination and proteasome degradation of Bcl-2 in the human acute lymphoblastic leukemia cell line, Reh. *Int. J. Mol. Med.*, **5**, 165-171 (2000).
6. Bai, X.C., et al., Phospholipase C γ_1 is required for survival in heat stress: involvement of protein kinase C-dependent Bcl-2 phosphorylation. *J. Biochem.*, **131**, 207-212 (2002).
7. Jarvis, W.D., et al., Evidence for involvement of mitogen-activated protein kinase, rather than stress-activated protein kinase, in potentiation of 1-beta-D-arabinofuranosylcytosine-induced apoptosis by interruption of protein kinase C signaling. *Mol. Pharmacol.*, **54**, 844-856 (1998).
8. Wall, N.R., et al., Modulation of cIAP-1 by novel antitubulin agents when combined with bryostatin 1 results in increased apoptosis in the human early pre-B acute lymphoblastic leukemia cell line Reh. *Biochem. Biophys. Res. Commun.*, **9**, 76-80 (1999).

AH/LY/MAM 1/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.