

Product Information

XIAP

Human, Recombinant
Expressed in *Escherichia coli*

Product Number **X 3378**

Product Description

XIAP is produced from a DNA sequence corresponding to amino acids 1-497 of human XIAP fused with the amino acid sequence (MATVIDHHHHHSSNG) at the amino terminus. By immunoblotting, recombinant human full length XIAP migrates as a 58 kDa protein under reducing and non-reducing conditions.

XIAP (X-linked inhibitor of apoptosis protein) is a member of the IAPs (inhibitor of apoptosis proteins) that function in cell death pathways to inhibit programmed cell death.¹ XIAP functions downstream of mitochondrial events to inhibit apoptosis. The IAPs share one to three copies of an approximately 70 amino acid sequence motif, BIR (baculovirus IAP repeat).² The BIR domains promote protein-protein interactions with caspases as well as with members of the TRAF family of signal molecules.³ The BIR2 domain of XIAP inhibits caspase-3 and caspase-7⁴⁻⁶ and the BIR3 domain inhibits caspase-9.^{7,8} The ability of XIAP to inhibit caspases is prevented by SMAC-Diablo through binding to XIAP-BIR2 and BIR3.⁹ IAPs are widely but differentially expressed in human cancers and leukemias.¹⁰

Reagent

Recombinant Human XIAP is supplied at approximately 0.224 mg/ml from a 0.2 μ m filtered solution in 25 mM HEPES (pH 8.0), 0.1 M sodium chloride, and 1 mM dithiothreitol (DTT).

Storage/Stability

Store at -20 °C. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. Do not store in a frost-free freezer.

Product Profile

Inhibition of the fluorogenic caspase substrate, DEVD-afc (N-acetyl-Asp-Glu-Val-Asp-7-amido-4-trifluoromethylcoumarin, Sigma Product No. A 0466) cleavage activities in cell extracts (2×10^6 cells) activated by the addition of cytochrome c and dATP, is 100 pM for recombinant human caspase-3 and by 350 pM for recombinant human caspase-7

The typical amount of XIAP required to decrease the rate of DEVD-afc cleavage by 50% in activated cell extracts is between 40-100 nM.

The typical amount of XIAP for inhibition of 50% of recombinant human caspase-3 is between 10 and 200 nM.

The typical amount of XIAP for inhibition of 50% of recombinant human caspase-7 is between 2 and 50 nM.

Note: In order to obtain the best results in various techniques and preparations, we recommend that each laboratory should determine optimal working dilutions.

Purity: >95% as determined by SDS-PAGE, visualized by silver stain.

References

1. Wagenknecht, B., et al., *Cell Death Differ.*, **6**, 370-376 (1999).
2. Sun, C., et al., *Nature*, **401**, 818-822 (1999).
3. Roy, N., et al., *EMBO J.*, **16**, 6914-6925 (1997).
4. Deveraux, Q.L., et al., *EMBO J.*, **18**, 5242-5251 (1999).

5. Suzuki, Y., et al., J. Biol. Chem., **276**, 27058-27063 (2001).
6. Riedl, S.J., et al., Cell, **104**, 791-8000 (2001).
7. Ekert, P.G., et al., J. Cell Biol., **152**,483-490 (2001).
8. Deveraux, Q.L., et al., EMBRO. J., **17**, 2215-2223 (1998).
9. Verhagen, A.M., et al., Cell, **102**, 43-53 (2000).
10. Tamm, I., et al., Clin. Cancer Res., **6**, 1796-1803 (2000).

KAA 01/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.