

Milliflex[®] Quantum: A Fluorescence-based Platform for the **Rapid Detection of Contaminants in Filterable Products** Anne Baumstummler, Renaud Chollet, Hervé Meder, Sophie Rouillon, Frédéric Olivieri, Gaël Waiche and Sébastien Ribault

INTRODUCTION



Testing for microbial contamination in pharmaceutical production processes is necessary to monitor product quality. However, traditional microbiological methods are slow and require several days to obtain results. Therefore there is an increasing demand for rapid microbiological methods. Millipore's Milliflex Quantum system is a convenient and sensitive platform for the quantitative detection of contaminants in filterable samples. This rapid microbiological method is based on an universal enzymatic fluorescent staining of viable and culturable microorganisms. The fluorescent staining procedure is nondestructive, allowing downstream specific identification following a positive result. It was demonstrated that the Milliflex Quantum protocol is a fast and reliable alternative for the detection of microorganisms in water and complex matrices. The time to results were shorter than the time to results obtained with the compendial method.

PRINCIPLE OF DETECTION

fluorescence inside cells is an indicator of microbial metabolism The principle of the fluorescence detection is based on an enzymatic reaction. The fluorogenic substrate used is a non-fluorescent activity and membrane integrity. The dye is diluted in a staining viability marker that is cleaved by nonspecific intracellular buffer allowing cell membrane permeabilization and thus dye enzymes resulting in a fluorescent product. Accumulation of introduction into cells.

MATERIAL AND METHODS

(A) Sterile water artificially contaminated with ATCC or environmethod was performed in parallel to control results concordance mental germs or (B) inprocess nonsterile water samples from between both methods. Traditional plate count was done after a pharmaceutical plants were filtered over mixed cellulose ester 7- day incubation. Fluorescent microcolonies counts and visual membranes using the Milliflex PLUS pump system. After incubation counts after reincubation of stained membranes were compared on Milliflex cassettes prefilled with Tryptic Soy Agar (TSA) or R2A to results from compendial method. agar at 32.5 ± 2.5 °C, or Sabouraud Dextrose Agar (SDA) at 22.5 (C) Chinese Hamster Ovary (CHO) cells culture from 4 to 6.10^{6} ± 2.5 °C, membranes were transferred onto a cellulose pad soaked cells/mL was spiked with Bacillus cereus, Staphylococcus with the Milliflex Quantum staining solution and incubated for 30 epidermidis or Propionibacterium acnes and incubated at least minutes at 32.5 ± 2.5 °C. Fluorescent microcolonies were counted for 16 hours at 37 °C. The culture was then filtered over mixed using the Milliflex Quantum reader (specific LED system). After cellulose ester membranes and incubated on Milliflex cassettes detection, membranes were reincubated on media for traditional prefilled with TSA at 37 °C. After the incubation period, memplate count (viability assay) and contaminants identification branes were stained following the procedure described above. using the MicroSEQ[®] platform (Applied Biosystems). Compendial



RESULTS



Microorganisms		Incubation Condition	Incubation time required for Milliflex Quantum detection
Escherichia coli	ATCC 8739	TSA 32.5 ± 2.5°C	8 hours
Pseudomonas aeruginosa	ATCC 9027		12 hours
Bacillus subtilis	ATCC 6633		8 hours
Staphylococcus aureus	ATCC 6538		12 hours
Candida albicans	ATCC 10231	SDA	24 hours
Aspergillus brasiliensis	ATCC 16404	22.5 ± 2.5°C	30 hours
Caulobacter spp.	Environmental strains	R2A 32.5 ± 2.5°C	30 hours
Micrococcus Iylae			16 hours
Ralstonia spp.			24 hours

B Inprocess Nonsterile Water Samples from Pharmaceutical Plants





Milliflex Quantum

Detection

Fluorescent Microcolonies 40 hours

Visual Plate after Reincubation



Contaminants Identification: Rhodococcus spp.

C Detection of Contaminants in Mammalian Cell Culture

Visible Colonies

7 days

Incubation time (on TSA at 37°C) required for the detection of different bacteria with the Milliflex Quantum system



8 hours

S. epidermidis 9 hours



P. acnes 48 hours

CONCLUSION

Spiked or environmental germs were efficiently and rapidly detected with the Milliflex Quantum protocol with a 3 to 5 times reduction of the incubation time compared to the compendial method. As the method is nondestructive, each fluorescent microcolony detected will continue to grow to yield visible colonies allowing the identification of the contaminants using any available identification method. The study demonstrated that the Milliflex Quantum platform could easily replace the traditional microbiological method with up to 5 days saved for slow growing strains. This method was also efficiently applied to the detection of contaminants in mammalian cell culture showing that the Milliflex Quantum system can be easily adapted to complex matrices. Detecting contamination early allows production to implement corrective action sooner.





Millipore and Milliflex are registered trademarks of Millipore Corporation. The M mark is a trademark of Millipore Corporation. MicroSEQ is a registere trademark of Applied Biosystems, LLC. Lit. No. PS2734EN00

© 2010 Millipore Corporation, Billerica, MA 01821 U.S.A. All rights reserved.