

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

MERCK

CONCENTRATE MORE, WORRY LESS

with **ULTRA High Recoveries**

Achieve the fastest and most efficient concentration for even the most sensitive downstream applications, thanks to these unique features of AMICON® Ultra filters:

- **High sample recovery:** 90% recovery from even the most dilute samples thanks to the low analyte binding of our Ultracel® membranes
- **Reliable:** a true dead stop provides consistent recovery volumes without spinning to dryness
- **Less clogging, faster filtering:** the dual vertical membrane design is aligned with the filtrate rather than perpendicular to it
- **Quick processing time:** 10–30 minute spin times
- **Optimal sample retrieval:** the reverse spin feature allows for maximum recovery, especially important with dilute samples

Ultra available!

To learn more, please visit:
SigmaAldrich.com/AMICON



AMICON® Ultra offer* 25% off

Product Description	Cat No.
AMICON® ULTRA 0.5 mL 3K 8PK	UFC500308
AMICON® ULTRA 0.5 mL 3K 24PK	UFC500324
AMICON® ULTRA 0.5 mL 3K 96PK	UFC500396
AMICON® ULTRA 0.5 mL 10K 8PK	UFC501008
AMICON® ULTRA 0.5 mL 10K 24PK	UFC501024
AMICON® ULTRA 0.5 mL 10K 96PK	UFC501096
AMICON® ULTRA 0.5 mL 30K 8PK	UFC503008
AMICON® ULTRA 0.5 mL 30K 24PK	UFC503024
AMICON® ULTRA 0.5 mL 30K 96PK	UFC503096
AMICON® ULTRA 0.5 mL 50K 8PK	UFC505008
AMICON® ULTRA 0.5 mL 50K 24PK	UFC505024
AMICON® ULTRA 0.5 mL 50K 96PK	UFC505096
AMICON® ULTRA 0.5 mL 100K 8PK	UFC510008
AMICON® ULTRA 0.5 mL 100K 24PK	UFC510024
AMICON® ULTRA 0.5 mL 100K 96PK	UFC510096
AMICON® ULTRA 2 mL 3K 24PK	UFC200324
AMICON® ULTRA 2 mL 10K 24PK	UFC201024
AMICON® ULTRA 2 mL 30K 24PK	UFC203024
AMICON® ULTRA 2 mL 50K 24PK	UFC205024
AMICON® ULTRA 2 mL 100K 24PK	UFC210024
AMICON® ULTRA-4 4 mL 3K NMWL 8PK	UFC800308
AMICON® ULTRA-4 4 mL 3K NMWL 24PK	UFC800324
AMICON® ULTRA-4 4 mL 3K NMWL 96PK	UFC800396
AMICON® CENTRIFUGAL 4 mL 10K NMWL 8pk	UFC801008
AMICON® CENTRIFUGAL 4 mL 10K NMWL 24pk	UFC801024

Product Description	Cat No.
AMICON® CENTRIFUGAL 4 mL 10K NMWL 96pk	UFC801096
AMICON® CENTRIFUGAL 4 mL 30K NMWL 8pk	UFC803008
AMICON® CENTRIFUGAL 4 mL 30K NMWL 24pk	UFC803024
AMICON® CENTRIFUGAL 4 mL 30K NMWL 96pk	UFC803096
AMICON® ULTRA-4 50k NMWL 8pk	UFC805008
AMICON® ULTRA-4 50k NMWL 24pk	UFC805024
AMICON® ULTRA-4 50k NMWL 96pk	UFC805096
AMICON® ULTRA-4 100k NMWL 8pk	UFC810008
AMICON® CENTRIFUGAL 4 mL 100K NMWL 24pk	UFC810024
AMICON® CENTRIFUGAL 4 mL 100K NMWL 96pk	UFC810096
AMICON® ULTRA 15 mL 3K NMWL 8PK	UFC900308
AMICON® ULTRA 15 mL 3K NMWL 24PK	UFC900324
AMICON® ULTRA 15 mL 3K NMWL 96PK	UFC900396
AMICON® ULTRA 15 mL 10K NMWL 8PK	UFC901008
AMICON® ULTRA 15 mL 10K NMWL 24PK	UFC901024
AMICON® ULTRA 15 mL 10K NMWL 96PK	UFC901096
AMICON® ULTRA 15 mL 30K NMWL 8PK	UFC903008
AMICON® ULTRA 15 mL 30K NMWL 24PK	UFC903024
AMICON® ULTRA 15 mL 30K NMWL 96PK	UFC903096
AMICON® ULTRA 15 mL 50K 8PK	UFC905008
AMICON® ULTRA 15 mL 50K 24PK	UFC905024
AMICON® ULTRA 15 mL 50K 96PK	UFC905096
AMICON® ULTRA 15 mL 100K 8pk	UFC910008
AMICON® ULTRA 15 mL 100K 24pk	UFC910024
AMICON® ULTRA 15 mL 100K 96pk	UFC910096

Use promo code **IBW** to have access to special pricing

***Limited Time Offer! Hurry Up!**

WANT TO STAY UP-TO-DATE
with all upcoming seasonal offers?

Subscribe to our email reminder and don't miss one of the upcoming attractive promotion and prizes.

**SigmaAldrich.com/
SeasonalSpecials**



To place an order or receive technical assistance

Order/Customer Service: SigmaAldrich.com/order
Technical Service: SigmaAldrich.com/techservice

SigmaAldrich.com

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt, Germany

MerckMillipore.com

We have built a unique collection of Life Science brands with unrivalled experience in supporting your scientific advancements.

Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®

© 2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, the vibrant M, Millipore, Sigma-Aldrich, and AMICON® are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

MK_FL13607EN Ver. 1.0 55575 05/2024

