



CHICKEN ANTI-MONOCARBOXYLATE TRANSPORTER 2 (MCT2) POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB1287
LOT NUMBER:	
QUANTITY:	50 µg
CONCENTRATION:	xx mg/mL
SPECIFICITY:	Monocarboxylate (lactate) transporter 2
IMMUNOGEN:	Synthetic peptide with the amino acid sequence NTHNPPSDRDKESSI corresponding to the carboxyl-terminal of rat MCT2.
APPLICATIONS:	Western Blot: 1:10,000 Immunohistochemistry on tissue fixed in neutral buffered formaldehyde: 1:300 Works on paraffin embedded tissue sections. Immunocytochemistry Immunoelectron microscopy Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITY:	Rat
FORMAT:	Affinity purified chicken IgY.
PRESENTATION:	Liquid in PBS and 35% glycerol. Contains no preservative.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.
RELATED REFERENCES:	<i>Cell</i> (1994) 76 :865-873. <i>J. Biol. Chem.</i> (1995) 270 :1843-1849. <i>Am. J. Physiol.</i> (1993) 264 :C761-C782.
REFERENCE:	<i>Glia</i> (1998) 22 :272-281. <i>J. Biol. Chem.</i> (1999) 274 :28420-28426.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2010: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com