

QUALITY CONTROL RANGES

Milliplex Mouse Cytokine Magnetic Bead Panel 2 Catalog# MCYP2MAG-73K

Control Catalog # MXM6073

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
*EPO	Control 1	MCYP2-105		pg/mL
	Control 2	MCYP2-205		pg/mL
*Exodus-2	Control 1	MCYP2-105		pg/mL
	Control 2	MCYP2-205		pg/mL
MCP-5	Control 1	MCYP2-105	135 – 281	pg/mL
	Control 2	MCYP2-205	692 - 1437	pg/mL
MIP-3 β	Control 1	MCYP2-105	146 - 304	pg/mL
	Control 2	MCYP2-205	709 - 1473	pg/mL
TARC	Control 1	MCYP2-105	141 - 293	pg/mL
	Control 2	MCYP2-205	705 - 1465	pg/mL
MIP-3a	Control 1	MCYP2-105	167 - 347	pg/mL
	Control 2	MCYP2-205	689 - 1431	pg/mL
*IL-16	Control 1	MCYP2-105		pg/mL
	Control 2	MCYP2-205		pg/mL
*Fractalkine	Control 1	MCYP2-105		pg/mL
	Control 2	MCYP2-205		pg/mL
IL-21	Control 1	MCYP2-105	119 - 246	pg/mL
	Control 2	MCYP2-205	649 - 1364	pg/mL
IL-22	Control 1	MCYP2-105	147 - 306	pg/mL
	Control 2	MCYP2-205	711 - 1477	pg/mL
IL-25	Control 1	MCYP2-105	1480 - 3074	pg/mL
	Control 2	MCYP2-205	7150 - 14851	pg/mL
IL-28B	Control 1	MCYP2-105	152 - 315	pg/mL
	Control 2	MCYP2-205	758 - 1573	pg/mL

EMD Millipore • (636) 441-8400 • (636) 441-8050 (FAX) • www.millipore.com

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

***QC ranges were not set for these analytes.**