

QUALITY CONTROL RANGES
MILLIPLEX® MAP Human Pituitary Magnetic Bead Panel 1 Catalog# HPTP1MAG-66K



MILLIPLEX® MAP Non-Human Primate Pituitary Magnetic Bead Panel 1 Catalog# NHPPT1MG-46K

Control Catalog # PIT-6046
Lot # PIT-109 & PIT-209

<u>Hormone</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>	<u>Hormone</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
ACTH	Control 1	44 - 91	pg/mL	HGH	Control 1	38 - 79	pg/mL
	Control 2	175 - 364	pg/mL		Control 2	157 - 325	pg/mL
AGRP	Control 1	32 - 67	pg/mL	HLH	Control 1	0.62 - 1.3	mIU/mL
	Control 2	134 - 279	pg/mL		Control 2	2.5 - 5.3	mIU/mL
CNTF	Control 1	1552 - 3224	pg/mL	HTSH	Control 1	0.53 - 1.1	μIU/mL
	Control 2	6817 - 14159	pg/mL		Control 2	2.2 - 4.5	μIU/mL
RFSH	Control 1	0.32 - 0.66	mIU/mL				
	Control 2	1.3 - 2.8	mIU/mL				

Note: The Quality Control Ranges are generated with using serum matrix provided in the kit. Quality Control values in culture media are not tested.

PIT-6046: PIT-109.209 (15 March 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.