

## Product Information

### Aldrich 744077 COware

Screwable two-chamber glass system

### Warning!

Sigma-Aldrich cannot guarantee glass equipment from breakage under pressure because of varying conditions.

- Glass equipment should always be examined for damages to its surface, which may weaken its strength.
- One must abide to all laboratory safety procedures and always work behind a shield when working with glass equipment under pressure.
- Materials which solidify upon standing thus causing stress on the glass equipment should be avoided.
- The equipment should not be operated above 60 psi/4 bar.

