



CHICKEN ANTI-NEURON SPECIFIC ENOLASE [NSE] AFFINITY PURIFIED POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB9698	QUANTITY:	500 μ L
LOT NUMBER:		CONCENTRATION:	XX mg/mL
SPECIFICITY:	Recognizes Neuron Specific Enolase (NSE).		
IMMUNOGEN:	Synthetic peptides from human/mouse/rat Neuron Specific Enolase.		
APPLICATIONS:	Immunohistochemistry: 1:1000-1:5000 using 2% paraformaldehyde fixed tissues. Immunocytochemistry: 1:1000-1:5000 using 2% paraformaldehyde fixed cells. Western blotting: On mouse brain extracts, blot shows a single band at 47 kDa. Optimal working dilutions must be determined by the end user.		
SPECIES REACTIVITY:	Human, mouse, and rat. Other species have not been tested.		
FORMAT:	Affinity purified chicken IgY. Two different antibodies were combined to make this product.		
PRESENTATION:	Liquid in PBS, pH 7.2 with 0.02% sodium azide.		
STORAGE/HANDLING:	Maintain at 2-8°C in undiluted aliquots for up to 6 months after date of receipt.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Also Available - Chicken IgY Negative Control Antibody (Cat. No. PP25),
Chicken IgY Blocking Solution (Cat. No. 2180) and Chicken IgY
Precipitating
Agarose (Cat. No. 2185)

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com