



QC RANGE SHEET
MILLIPLEX[®] MAP Human Neuropeptide Magnetic Bead Panel
Kit Catalog # HNP-MAG-35K

QC Catalog #:
Lot# HNP-111 and HNP-211

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
α -MSH	Control 1	231 – 480	pg/mL
	Control 2	921 – 1913	pg/mL
β -Endorphin	Control 1	574 – 1192	pg/mL
	Control 2	2370 – 4921	pg/mL
Neurotensin	Control 1	319 – 662	pg/mL
	Control 2	1293 – 2686	pg/mL
Oxytocin	Control 1	474 – 984	pg/mL
	Control 2	1827 – 3794	pg/mL
Substance P	Control 1	56 – 115	pg/mL
	Control 2	232 – 483	pg/mL

HNP-6035

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.

HNP-6035: HNP-111.211 (04 June 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA