

Product Information

Anti-Nucleolin

produced in rabbit, affinity isolated antibody

Catalog Number **N2662**

Product Description

Anti-Nucleolin is produced in rabbit using as immunogen a synthetic peptide corresponding to amino acid residues 2-17 of human nucleolin with C-terminal added cysteine, conjugated to KLH. The corresponding sequence differs by 4 amino acids in rat and mouse nucleolin and by one amino acid in the rat homolog protein NRP (nucleolin related protein). The antibody is affinity-purified using the immunizing peptide immobilized on agarose.

Anti-Nucleolin recognizes human, mouse, and rat nucleolin. Applications include immunoblotting (105 kDa), immunoprecipitation, and immunofluorescence. Detection of the nucleolin band by immunoblotting is specifically inhibited with the immunizing peptide.

Nucleolin, also known as C23, is the major nucleolar protein of exponentially growing eukaryotic cells. It is a ubiquitous conserved multi-domain phosphoprotein found in all vertebrate species.¹⁻³ Nucleolin contains a specific bipartite nuclear localization signal sequence and three different structural domains. The N-terminal domain comprises of highly acidic regions separated by basic sequences and contains multiple phosphorylation sites. The central domain contains four RNA-binding domains (RBD or RRM) and the C-terminal domain (GAR or RGG) is rich in glycine, arginine, and phenylalanine residues. It was suggested that the N- and C-terminal domains are involved in protein-protein interaction whereas the central domain is involved in specific interactions with nucleic acids.³

Nucleolin is a multifunctional protein involved in most steps of ribosome biogenesis. Nucleolin has been implicated in chromatin structure, rDNA transcription, rRNA maturation, ribosome assembly, cytokinesis, nucleogenesis, cell proliferation and growth, and nucleo-cytoplasmic transport. Nucleolin also exhibits autodegradation, DNA and RNA helicase activities, and DNA-dependent ATPase activity. Nucleolin activities may be regulated by proteolysis, methylation, ADP-ribosylation, and phosphorylation by several kinases, including casein kinase II (CK2), p34^{cdc2}, and protein kinase C- ζ .^{3,4}

Although nucleolin is mainly localized in the nucleus, several reports indicate its presence at the cell surface of cultured tumor cells and endothelial cells of angiogenic vessels *in vivo*.⁵⁻⁸ Nucleolin might function as a cell surface receptor for lipoproteins, viruses, extracellular matrix, growth factors, and other molecules.³⁻⁴ It has been found that nucleolin interacts with telomerase and that this interaction is critical for the nucleolar localization of telomerase.⁹

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

Antibody Concentration: 1.0-1.5 mg/mL

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing, or storage in frost-free freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilutions should be discarded if not used within 12 hours.

Product Profile

Immunoblotting: a working antibody concentration of 0.5-1.0 $\mu\text{g}/\text{mL}$ is recommended using a whole extract of mouse NIH-3T3 cells and rat PC12 cells.

Immunoprecipitation: 1-2 μg of the antibody immunoprecipitates nucleolin from 250 μg RIPA lysate of human HeLa cells.

Indirect immunofluorescence: a working antibody concentration of 10-20 $\mu\text{g}/\text{mL}$ is recommended using human HeLa cells.

Note: In order to obtain the best results using various techniques and preparations, we recommend determining the optimal working dilutions by titration.

References

1. Lapeyre, B., et al., *Proc. Natl. Acad. Sci. USA*, **84**, 1472-1476 (1987).
2. Srivastava, M., et al., *FEBS Lett*, **250**, 99-105 (1989).
3. Ginisty, H., et al., *J. Cell Sci.*, **112**, 761-772 (1999).
4. Srivastava, M., and Pollard, H.B., *FASEB J.*, **13**, 1911-1922 (1999).
5. Deng, J.S., et al., *Mol. Biol. Rep.*, **23**, 191-195 (1996).
6. Said, E. A., et al., *J. Biol. Chem.*, **277**, 37492-37502 (2002).
7. Sinclair, J.F., and O'Brien, A.D., *J. Biol. Chem.*, **277**, 2876-2885 (2002).
8. Christian, S., et al., *J. Cell Biol.*, **163**, 871-878 (2003).
9. Khurts, S., et al., *J. Biol. Chem.*, **279**, 51508-51515 (2004).

DS,PHC 07/14-1