

Product Information

ANTI-MAP KINASE ACTIVATED PROTEIN KINASE-3 (MAPKAPK-3)

Developed in Sheep, Affinity Isolated Antibody

Product Number **M7681**

Product Description

Anti-MAP Kinase Activated Protein Kinase-3 (MAPKAPK-3) is developed in sheep using a synthetic peptide K-QAGSSSASQGCNNQ[G] that corresponds to the amino acids 368-382 of human MAPKAPK-3 conjugated to KLH/BSA as immunogen. Affinity isolated antibody is obtained by peptide immunoaffinity chromatography.

Anti-MAP Kinase Activated Protein Kinase-3 (MAPKAPK-3) reacts specifically with human MAPKAPK-3 (42 kD) by immunoprecipitation. The antibody does not react with MAPKAPK-2. The antibody does not interfere with kinase activity.

Anti-MAP Kinase Activated Protein Kinase-3 may be used for immunoprecipitation of active MAPKAPK-3 in cell lysates but is not recommended for immunoblotting.

The mitogen-activated protein kinase (MAPK) signal transduction cascade is characterized by a series of kinases, each of which subsequently phosphorylates and activates the next kinase in the cascade. For example, MAPK-activated protein kinase (MAPKAPK) is a serine/threonine kinase that is phosphorylated and activated by the MAP kinase homologue, stress-activated protein kinase-2 (SAPK-2). SAPK is activated in response to some form of cellular stress and can initiate adaptive responses, such as apoptosis. Following activation, MAPKAPK in turn phosphorylates substrates with the sequence H-X-R-X-X-S where H is any hydrophobic residue¹. Three different MAPKAPK's have been identified to date. MAPKAPK-1 differs in response to inhibitors and amino acid sequence². MAPKAPK-2 and -3 are both activated in response to cellular stress, interleukin-1 and tumor necrosis factor in KB and HeLa cells^{1,3}. Their amino acid sequence shows 75% homology to each other and both have similar structural features including a proline-rich N-terminus, a nuclear localization signal and conservation of key amino acid residues (Thr-222, Ser-272 and Thr-334)^{3,4}. It is thought that MAPKAPK-2 and -3 are likely to have either overlapping or identical substrates *in vivo*¹. Both MAPKAPK-2 and MAPKAPK-3 activity can

be inhibited with SB 203580 (Sigma Product No. S8307).

Reagents

The product is supplied as the affinity isolated antibody in 35mM Tris-glycine, pH 7.4, containing 30% glycerol.

Protein concentration is approximately 0.36 mg/ml by Bradford analysis.

Storage/Stability

Store at 0°C to -20°C. Aliquot to avoid repeated freezing and thawing. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

Procedure

1. Add 10 µg of Anti-MAPKAP Kinase 3 to a microcentrifuge tube.
2. Add 100 µl of a 1:1 (v/v) Protein G-Agarose (Sigma Product No. P2294) slurry that has been washed with ice-cold PBS.
3. Add 200 µl of ice-cold PBS.
4. Incubate for 30 minutes to 1 hour at RT or overnight at 4°C in ice-cold PBS, rotating the sample so the antibody and protein G mix thoroughly.
5. Pellet at 14,000 rpm for 15 seconds.
6. Remove the supernatant and then wash the Protein G-agarose twice with ice-cold Buffer A (50 mM Tris, pH 7.5, 1 mM EDTA, 1 mM EGTA, 0.5 mM Na₃VO₄, 0.1% 2-mercaptoethanol, 1% Triton X-100, 5 mM sodium pyrophosphate, 10 mM sodium glycerophosphate, 0.1 mM PMSF, 1 µg/ml aprotinin, 1 µg/ml leupeptin and 50 mM NaF). Note: Add PMSF fresh just before using buffer.
7. Resuspend the pellet of washed beads in 100 µl of ice-cold Buffer A.
8. Add sample containing antigen to beads.^a
9. Add 200 µl of ice-cold PBS
10. Incubate for 2 hours at 4°C, rotating the sample to mix thoroughly.

11. Wash the antigen/antibody/protein G-agarose complex with 500 μ l of ice-cold Buffer A containing 0.5 M NaCl. Repeat the wash.
12. After the second wash, remove the supernatant and resuspend the pellet in 20 to 50 μ l of Laemmli sample buffer (2% SDS, 10% glycerol, 0.05 M Tris, pH 6.8, containing bromphenol blue as the dye marker). Heat for 5 minutes at 100°C.
13. Microcentrifuge for 5 seconds. Transfer the supernatant to a fresh tube.
14. If a non-reducing gel is to be run, the sample can be loaded directly. If a reducing gel is to be run, add 5% 2-mercaptoethanol, incubate 1 hour 37°C then load sample onto gel.
15. After running the gel, analyze by protein staining, immunoblotting, or autoradiography (if the sample was radiolabeled).

^a In order to obtain the best results, we recommend trying several amounts of sample to a given amount of beads in order to determine the optimal condition for immunoprecipitation.

Product Profile

Recommended use: 5 μ g of Anti-MAP Kinase Activated Protein Kinase-3 will immunoprecipitate active MAPKAPK-3 from a cell lysate of HeLa cells stimulated with 0.5 mM sodium meta-arsenite.

Note: In order to obtain best results and assay sensitivity in different techniques and preparations we recommend determining optimal working dilutions by titration test.

References

1. Clifton, A.D., Young, P.R. and Cohen, P., *FEBS Lett.*, **392**, 209 (1996).
2. Stokoe, D., Campbell, D.G., Nakielny, S., Hidaka, H., Leever, S.J., Marshall, C. and Cohen, P. *EMBO J.*, **11**, 3985 (1992).
3. McLaughlin, M.M., Kumar, S., McDonnell, P.C., Van Horn, S., Lee, J.C., Livi, G.P. and Young, P.R., *J. Biol. Chem.*, **271**, 8488 (1996).
4. Sithanandam, G., Latif, F., Duh, F.M., Bernal, R., Smola, U., Li, H., Kuzmin, I., Wixler, V., Geil, L. and Shrestha, S., *Molec. Cell. Biol.*, **16**, 868 (1996).

coh 3/99

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.