



QUALITY CONTROL RANGES
MILLIPLEX[®] MAP Human CVD Panel 3 (Acute Phase) Magnetic Bead Panel
Catalog # HCVD3MAG-67K
Control Catalog # HCVMG3-6067-3 Lot # HCVM3-111 and HCVM3-211

<u>Hormone</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
A2M	Control 1	6.3 - 13	ng/mL
	Control 2	66 - 138	ng/mL
CRP	Control 1	0.14 - 0.29	ng/mL
	Control 2	1.4 - 3.0	ng/mL
Fetuin A	Control 1	3.3 - 6.9	ng/mL
	Control 2	33 - 69	ng/mL
AGP	Control 1	1.2 - 2.5	ng/mL
	Control 2	13 - 26	ng/mL
Fibrinogen	Control 1	0.32 - 0.66	ng/mL
	Control 2	3.1 - 6.4	ng/mL
L-Selectin	Control 1	0.05 - 0.10	ng/mL
	Control 2	0.47 - 1.0	ng/mL
SAP	Control 1	0.05 - 0.10	ng/mL
	Control 2	0.46 - 1.0	ng/mL
Haptoglobin	Control 1	0.77 - 1.6	ng/mL
	Control 2	7.5 - 17	ng/mL
PF4/CXCL4	Control 1	0.59 - 1.2	ng/mL
	Control 2	5.7 - 12	ng/mL
Adipsin	Control 1	1.1 - 2.4	ng/mL
	Control 2	12 - 26	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format.

HCVMG3-6067-3: HCVM3-111.211 (08 November 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA